RNA Isolation Protocol in Bixa orellana V.1

DOI

<https://dx.doi.org/10.17504/protocols.io.vy9e7z6>

Victor Carballo-Uicab¹, Yair Cárdenas-Conejo², Alba Vallejo-Cardona³, Margarita Aguilar-Espinosa¹,
Jacobo Rodríguez-Campos⁴, Hugo Serrano Posada², José Iberto Narváez Zapata⁵, Renata Rivera-Madrid¹

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Victor Carballo

cicy

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

DOI: <https://dx.doi.org/10.17504/protocols.io.vy9e7z6>

Protocol Citation: Victor Carballo-Uicab, Yair Cárdenas-Conejo, Alba Vallejo-Cardona, Margarita Aguilar-Espinosa, Jacobo Rodríguez-Campos, Hugo Serrano Posada, José Iberto Narváez Zapata, Renata Rivera-Madrid . RNA Isolation Protocol in Bixa orellana. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.vy9e7z6>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Created: November 29, 2018

Last Modified: November 29, 2018

Protocol Integer ID: 18177

Keywords: bixin, carotenoids, pigment, commercial rna extraction protocol, length cdna from bixa orellana, rna extraction, bixa orellana, long distance reverse transcript pcr protocol, rna, large quantities of pigment, length cdna

Abstract

A reliable protocol is described for isolation of large full-length cDNA from *Bixa orellana* mature tissues containing large quantities of pigments, phenols, and polysaccharides. This protocol involves the optimization of a commercial RNA extraction protocol in combination with a long distance reverse transcript PCR protocol.

Guidelines

Rodriguez-vila et al. - 2009 - Full-le...

Materials

liquid nitrogen

Pure-Link Micro-to-Midi kit was used to extract RNA (Cat. No. 18068-018 Invitrogen®)

DNase (DNase I Amp Grade, Deoxyribonuclease I Amplification Grade, Invitrogen, Cat. No 18068-015)

SuperScript™ III Reverse Transcriptase (Invitrogen, Cat. No. 18080-093)

- 1 Total RNA from *Bixa orellana* L. seeds was isolated from each of the five developmental stages in triplicate from three biological replication of N4P and P13W. Seeds ( 50 mg 50 mg) were ground in liquid nitrogen using a mortar and pestle.
- 2 A Pure-Link Micro-to-Midi kit was used to extract RNA (Cat. No. 18068-018 Invitrogen®), which was subsequently treated with DNase (DNase I Amp Grade, Deoxyribonuclease I Amplification Grade, Invitrogen, Cat. No 18068-015) and stored at  -80 °C until analysis. The RNA concentration was measured by spectrophotometric analysis (NanoDrop® ND-1000), and its quality was verified by agarose gel electrophoresis at  1.2 Mass Percent 1.2 % and ethidium bromide staining. Reverse transcription was performed with SuperScript™ III Reverse Transcriptase (Invitrogen, Cat. No. 18080-093) using  100 ng of total RNA from each tissue.